

**Product Name :**  
GYRATORY Sieve Shaker

**Product Code :**  
NLAB-TESTINGLABM109072



**Description :**

GYRATORY Sieve Shaker

**Technical Specification :**

Carries upto 7 sieves of 200mm or 300mm, 450mm diameter. The shaker is driven by a  $\frac{1}{4}$  H.P. Motor through a reduction gear immersed in oil. The sieve table does not rotate but is inclined from the vertical axis and the direction on inclination charges progressively in clockwise direction. If the stop pin below the table is removed, the shaker can have a rotary motion. In addition to this gyratory motion of the table, there is an upward and downward movement ensuring that each square cm of the sieve is utilized. A pair of rods and a holder are supplied. The holder can be fixed on the top of the upper most sieve, and thus the sieve set is firmly held. Suitable for operation from 230 V, 50 Hz Single Phase A.C. Supply.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

