

Product Name :
Gyroscope

Product Code :
NLAB-TECHNICALAB35008



Description :

Gyroscope

Technical Specification :

Bench-top, to demonstrate free and forced precession. High speed flywheel (up to 7500 rpm) guarded and mounted in a cradle fitted with ball bearings, external moment weight with scale , main and precession motors, control gear and tachometers mounted in apparatus base. Power supply 230V/50Hz.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003