

**Product Name :**  
Haemocytometer Disposable

**Product Code :**  
NLAB-EDUSCHOOLEQ173019

**Description :**

Haemocytometer Disposable

**Technical Specification :**

Features and benefits:

Bright grid pattern

2 tests per slide

No need to coverslip

Accurate, reliable and repeatable

No more contact to hazardous materials

Sturdy and strong

Quartz grade optical plastic

No need to wash for reuse

Main applications:

Microbiology, bacterial and fungal spores

Cell culture, cell concentration measurement and cell viability counts

Blood analysis

Chamber volume : 10  $\mu$ L

Chamber depth: 100  $\mu$ m

Dimensions: 25 x 75 x 1.6 mm (w x l x t)

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003