

**Product Name :**  
Half Body CPR Training Model Male

**Product Code :**  
NLAB-MEDICALMODELAB11027



**Description :**

Half Body CPR Training Model Male

**Technical Specification :**

Simulate standard open airway

External breast compression, intensity display of correct (at least 5cm):

The yellow indicator light display of correct compression position

Buzzing sound of correct compression intensity, the green indicator light display

Alarming sound of wrong compression intensity, the red indicator light display

Artificial respiration (inhalation): the volume of inhalation is observed through observing the wave of breast;

Operation methods: exercise operation

Power supply:- Battery

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003