

| | |
|---|---|
| Product Name : Halogen Spot Light | Product Code : NLABMEDICAL-M3220014 |
| Description : Halogen Spot Light India, Halogen Spot Light manufacturer India, suppliers, Medical Halogen Spot Light India, Hospital Halogen Spot Light India | |
| Technical Specification : Halogen spotlight is equally useful in operation theater OPD department and clinics with 12V/50W halogen blub & 12V/50W variable transformer with low & high intensity control Floor stand model & also required of table clamp model | |

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003