

Product Name :
Hand Held Garmin GPS Units

Product Code :
NLAB-SCIENTIFICIND12003



Description :

Hand Held Garmin GPS Units

Technical Specification :

Used for preliminary surveys & mobile data collection and tracking of geo-references for student education, With precision of at least 5m, Physical and Performance

Dimensions: 2.3"x4.5"x1.4"

Display size: 1.53"x2.55"; 3" diagonal

Resolution: 240 x 400 pixels

Display type: transfective color TFT touchscreen

Battery: 2 AA (not included); NiMH or Lithium recommended

Waterproof: yes

Floats: no

High sensitivity receiver: yes

Interface: high-speed USB and NMEA 0183 compatible

Maps and Memory

Base map: yes

Preloaded maps: yes (topographic)

Ability to add maps: yes

Built-in memory: 850 MB

Accepts data cards: micro SD card (not included)

Features

Automatic routing: yes

Electronic compass: yes

Picture view: yes

Garmin Connect compatible: yes

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003