

| | |
|---|---|
| Product Name : Hand Held Pulse Oximeter | Product Code : NLABMEDICAL-M3410003 |
|---|---|

| |
|--|
| |
|--|

Description :

Hand Held Pulse Oximeter

Technical Specification :

Digital Technology

Alarm with freely set limit

Accurate measurement during low perfusion

Anti-Movement, Stable Performance

Built-in rechargeable battery and seat-type charger, Power-Saved Function

Built-in memory card for data memory, connection with computer for downloading data

Bright, Easy-to-Read LED Numerals

Light weight, Portable 2 parameters: SPO2 Pulse Rate

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003