

Product Name : Hand Held Pulse Oximeter	Product Code : NLABMEDICAL-M3410003
Description :	
Hand Held Pulse Oximeter	
Technical Specification :	
Digital Technology	
Alarm with freely set limit	
Accurate measurement during low perfusion	
Anti-Movement, Stable Performance	
Built-in rechargeable battery and seat-type charger, Power-Saved Function	
Built-in memory card for data memory, connection with computer for downloading data	
Bright, Easy-to-Read LED Numerals	
Light weight, Portable 2 parameters: SPO2 Pulse Rate	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003