

**Product Name :**  
Hand Held Single Channel ECG

**Product Code :**  
NLABMEDICAL-M3380001



**Description :**

Hand Held Single Channel ECG

**Technical Specification :**

Input circuit: Floating.

Calibrating Voltage :1mV.

Input current: <50mA.

CMRR: >103dB.

A/D: 12bit.

Frequency Response: 0.05Hz--150Hz.

Filter: AC 50/60Hz EMG: 35-45Hz.

Skin Voltage Tolerance: +300mV.

Sensitivity Setting 5, 10, 20mm/mV, conversion deviation 25%.

Noisy Level: ?15µVp-p.

Safety standard: IEC class I, type CF.

Lead selection: Automatic and Manual.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003