

Product Name : Hand Opearted Pipe Bender	Product Code : NLAB-MACHINELABWK25004
--	---



Description :
Hand Opearted Pipe Bender
Technical Specification :
Hand Opearted Pipe Bender .
 SPECIFICATIONS Working Round Plate 510 x 22 mm Bending Capacity upto 25.4 mm x 8 swg. Working Height 1016mm Overall Dimension Length x Breadth x Height 406 mm x 406 mm x 1016 mm
Working Round Plate
Bending Capacity upto 25.4 mm x 8 swg.
Working Height
Overall Dimension Length 406 mm x Breadth 406 mm x Height 1016 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003