

Product Name :
Hand Operated Casagrande Equipment

Product Code :
NLAB-VOCATIONALTVET18002



Description :

Hand Operated Casagrande Equipment

Technical Specification :

Use:-

Measuring the Liquid limit.

Hand operated on the counter.

A brass cup suspended from a carriage designed to control its drop exactly 1cm onto a hard rubber base.

Adjustable height of drop by a nut at the rear of the cup.

Casagrande grooving tool.

Rubber feet isolate base from work surface

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com