

Product Name :
Hand Operated Hydraulic Press

Product Code :
NLAB-MACHINELABWK14007



Description :

Hand Operated Hydraulic Press

Technical Specification :

The Press is provided within a 2 stage hand pump. The low pressure pump is used for rapid advance of piston at low pressure. The high pressure Pump is provided with a special material lining and is used for slow advance at high pressure.

The winch provided for raising and lowering the table is smooth running. During operation the table is fastened by means of pins. The press has been all welded Standard accessories and planned off V Block for straightening.

Between Table and Bed

Length of

mm mm

360

360 100

360 125

360 150

360 180

360 190

360 200

360 200.

~~300~~ 210

~~300~~ 225

~~300~~ 240

~~305~~ 240

~~306~~ 250

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003