

Product Name :
Handheld Shearer

Product Code :
NLAB-MACHINELABWK34014



Description :

Handheld Shearer

Technical Specification :

Type: Handheld All steel construction

Length of blade: 300mm

Capacity for flat steel: 6mm

Capacity for Round steel: 13mm

With fixing hole at base.

Standard Welding Equipment Laboratory Set Handheld Shearer Manufacturer, Supplier and Handheld Shearer Exporter in India

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments