

**Product Name :**  
Handheld Spectrometer

**Product Code :**  
NLAB-SCHOOLABS21006



**Description :**

Handheld Spectrometer

**Technical Specification :**

Separation of light by triangular prism glass for clear results, resulting in light banding, different spectra. large clear bar Well packed in a wooden box with lock.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003