

Product Name :
Handheld Spectrometer

Product Code :
NLAB-SCHOOLABS21006



Description :

Handheld Spectrometer

Technical Specification :

Separation of light by triangular prism glass for clear results, resulting in light banding, different spectra. large clear bar Well packed in a wooden box with lock.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003