

**Product Name :**  
Hanger for Slotted Weight Steel Nickel Plated

**Product Code :**  
NLAB-EDUSCHOOLEQ136002



**Description :**

Hanger for Slotted Weight Steel Nickel Plated

**Technical Specification :**

Hanger for Slotted Weight 10 Gram Steel Nickel Plated.

Hanger for Slotted Weight 100 Gram Steel Nickel Plated.

Hanger for Slotted Weight 50 Gram Steel Nickel Plated.

Made from Nickel Plated Steel.

Polished for slotted weights and general lab use.

Designed to fit the following slotted weights:

80-7 Brass 1 Gram

80-8 Brass 2 Grams

80-1 Brass 5 Grams

80-2 Brass 10 Grams

80-3 Brass 20 Grams

80-4 Brass 50 Grams

80-5 Brass 100 Grams

80-6 Brass 200 Grams

80-9 Brass 500 Grams

75-51 Nickel Plated Steel 5 Grams

75-52 Nickel Plated Steel 10 Grams

75-53 Nickel Plated Steel 20 Grams

75-54 Nickel Plated Steel 50 Grams

75-55 Nickel Plated Steel 100 Grams

75-56 Nickel Plated Steel 200 Grams

75-59 Nickel Plated Steel 500 Grams

75-60 Nickel Plated Steel 1000 Grams

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

