

Product Name :
Head Protection Thermo Guard 100

Product Code :
NLAB-SAFETYLABP14008



Description :

Head Protection Thermo Guard 100

Technical Specification :

Adjustable chin strap

Size- 53cm-60cm

4 point textile suspension with slip fit or ratchet fit adjustment.

Proprietary injection-molded fiberglass compound

Features

8X stronger than HDPE helmets yet lightweight.

Aesthetic look with smooth finish – injection moulded solid colour resin for a permanent, no-chip, no-peel finish.

Environment friendly as the material is re-cyclable.

More rugged, to sustain extreme conditions.

Greater impact resistance and crack resistance with cross ribs more reliable and durable.

Designed for use in industries where there is a hazard from high levels of radiant heat.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003