

Product Name :
Head Protection UI 102

Product Code :
NLAB-SAFETYLABP14011



Description :

Head Protection UI 102

Technical Specification :

Peak: 23 mm

Size: 53 cm - 61 cm

Weight: 320 Grams

Colours: White (W) , Grey (Gr), Blue (B) ,Yellow (Y) ,Red (R) , Green (G),Orange (O), Purple (P)

Ventilation: Not available

Material:-

Head band: Textile

Shell: HDPE (High density poly ethylene)

Cradle: 8 point Textile Suspension

Sweat band: In-built

Storage:-

Can be stored and transported in their original cartons at ambient temperature (0°C to + 30°C)

Do not store in direct/high heat or sunlight as this may distort the shell

May be stored in the dark for up to 5 years

Lifetime:-

Excessive wear and tear can considerably reduce the lifespan of the product

Cleaning:-

Clean using warm soapy water and soft cloth

Do not use solvents or abrasives

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003