

**Product Name :**  
Head Protection UI 1211

**Product Code :**  
NLAB-SAFETYLABP14009



**Description :**

Head Protection UI 1211

**Technical Specification :**

Colours: White (W) , Grey (Gr), Blue (B) ,Yellow (Y) ,Red (R) , Green (G),Orange (O), Purple (P)

Peak: 48 mm

Ventilation: Not available

Accessories: Compatible with face accessories

Size: 53 cm - 61 cm

Weight: 320 Grams approx

Material:-

Shell: HDPE (High density poly ethylene)

Cradle: 6 point Plastic Suspension

Head band: Plastic

Sweat band: Not available

Storage:-

Can be stored and transported in their original cartons at ambient temperature (0°C to + 30°C)

May be stored in the dark for up to 5 years

Do not store in direct/high heat or sunlight as this may distort the shell

Lifetime:-

Excessive wear and tear can considerably reduce the lifespan of the product

Cleaning:-

Do not use solvents or abrasives

Clean using warm soapy water and soft cloth

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003