

Product Name :
Heat Insulation Tester

Product Code :
NLAB-TESTINGLABM93001



Description :

Heat Insulation Tester

Technical Specification :

Dimensions of sand bath : 400 x 400 mm (± 5 mm)

Temperature of sand bath : Adjustable up to 250°C

Diameter of copper disc : 15 ± 1 mm

Thickness of copper disc : 2.0 ± 0.1 mm

Diameter of steel balls : 5.0 ± 0.1 mm

Mass of steel balls : 4 kg (minimum)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

