

**Product Name :**  
Heat Resistance Tester Machines

**Product Code :**  
NLAB-TESTINGLABM107001



**Description :**

Heat Resistance Tester, Heat Resistance Tester Machines

**Technical Specification :**

Heat resistance tester. The heat resistance tester is ideal for assessing the effects of heat on outsole and upper materials. Heat is applied to the sample via a heated square pad manufactured in copper

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003