

Product Name :
Heat Resistance Tester Machines

Product Code :
NLAB-TESTINGLABM107001



Description :

Heat Resistance Tester, Heat Resistance Tester Machines

Technical Specification :

Heat resistance tester. The heat resistance tester is ideal for assessing the effects of heat on outsole and upper materials. Heat is applied to the sample via a heated square pad manufactured in copper

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003