

Product Name :
Heat Transfer Lab Equipments

Product Code :
NLAB-MECHANICALAB270001



Description :

Heat Transfer Lab Equipments

Technical Specification :

Current: 2 Amp

Material: Stainless Steel

Voltage: 220 V

Display: Digital

Current Type: DC

Vertical & Horizontal Condenser

Thermal Conductivity Of Insulating Slab

Unsteady State Heat Transfer Unit

Heat Transfer Through Lagged Pipe

Emissivity Measurement Apparatus

Heat Transfer Through Composite Walls

Finned Tube Heat Exchanger

Single Effect Evaporator

Thermal Conductivity of Liquids

Heat Pipe Demonstrator

Heat Transfer In Natural Convection

Thermal Conductivity Of Metal Rod

Dropwise Filmwise Condensation Apparatus

Heat Transfer In Forced Convection

Helical Coil Apparatus

Multi Effect Evaporator Double Effect Triple Effect

Stefan Boltzman Apparatus

Parallel Counter Flow Heat Exchanger

Shell & Tube Heat Exchanger

Open Pan Evaporator

Thermal Conductivity Of Insulating Powder

Heat Transfer From A Pin Fin

Heat Transfer In Agitated Vessel

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003