

**Product Name :**  
High Pressure Surgical Deluxe Autoclave

**Product Code :**  
NLAB-LABORATORYEM29003



**Description :**

High Pressure Surgical Deluxe Autoclave

**Technical Specification :**

High Pressure Surgical Deluxe Autoclave Ideal for all applications requiring, routine total destruction of all living microorganisms. Ideal for Clinics, Primary Health Centres, Nursing Homes, Hospitals, Doctors and Laboratories. Inner chamber & lid are made of thick stainless steel with pressure control switch, alarm & timer etc. Autoclave can be operated at any selected point 5 to 20 psi (pound per sq. inch).

Double Drum Size :

225mm x 300mm.

300mm x 375mm.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

