

Product Name :
High Pressure Vertical Autoclave Triple Walled Semi

Product Code :
NLAB-LABORATORYEM29008



Description :

High Pressure Vertical Autoclave Triple Walled Semi

Technical Specification :

Single panel feather touch LED microprocessor controller

Simultaneous display of temperature & time.

Timer starts automatically as soon as the set temperature is attained & stop the cycle automatically as soon the time is over.

Automatic priming : The air valve remains open upto 100°C for priming. As soon the unit achieves temp. of 100°C, the Solonoid Valve gets closed and starts generating pressure.

As soon as the cycle is complete, the steam is released automatically through solonoid valve.

Water may be finish & water level maintained automatically (Optional)

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com