

Product Name :
High Voltage Tester

Product Code :
NLAB-ELECTRICALAB160001



Description :

High Voltage Tester

Technical Specification :

High Voltage Tester to be used for a set of highervoltage up to 600kV.

The output should be single phase 50/60Hz .

High Voltage Tester

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003