

Product Name :
Hoek Triaxial

Product Code :
NLAB-ROCKTLAB170013



Description :

Hoek Triaxial

Technical Specification :

Hoek Cells have been designed to be used for triaxial testing of rock specimens.

Hoek

Cells comprise a steel body complete with two quick release self-sealing couplings, two steel end caps which are screwed to the cell body, 2 pieces of upper and 2 pieces of lower loading caps with spherical coupling and a rubber sealing sleeve to separate the specimen from the cell fluid.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003