

<b>Product Name :</b> Hoffman Voltameter	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730210
<b>Description :</b>	
Hoffman Voltameter	
<b>Technical Specification :</b>	
<p>The glass unit has two connected limbs, integral with reservoir tube and funnel shaped bulb, with a stopcock at the top of each side limb.</p> <p>The electrodes, available separately, are mounted in rubber stopper for insertion into the open ends of the Voltameter limbs.</p> <p>Pair of Carbon electrodes with brass terminals, fitted with rubber stoppers, approx 6mm dia x 100mm long in protection tubes.</p> <p>Metal stand with mounts, for giving support to voltmeter glass unit.</p> <p>Pair of platinum electrodes fused in glass tubes, fitted with rubber stoppers and protection tubes.</p> <p>Used on 12 V D.C. supply.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003