

**Product Name :**  
Hooked Weight Double Set

**Product Code :**  
NLAB-EDUSCHOOLEQ136001



**Description :**

Hooked Weight Double Set

**Technical Specification :**

- Double Hooked weight Set of 9, to ease in experiments when additional weights are required to be added.
- All weights have counter sunk recessed bottoms so that the weight rests flat.
- All weights have Brass Hooks.
- 1~1000g Steel, Black Painted
- 500g Steel, Black Painted
- 2~200 g Brass
- 1~100g Brass
- 1~50g Brass
- 2~20g Brass

- 1~10g Brass

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003