

**Product Name :**  
Horizontal Cylindrical Model

**Product Code :**  
NLAB-LABORATORYEM29004



**Description :**

Horizontal Cylindrical Model

**Technical Specification :**

Horizontal Cylindrical Model Inner chamber jacket and boiler is made of stainless steel. Sterilization cycle is controlled by three way valve or multi port valve stainless steel. It is fitted with automatic pressure control switch, indfoss piezostat, safety valve, pressure gauge, pressure release valve, water level indicator with cord & plug. Lid is made of mild steel lined with S.S. duly polished.

Size

Dia x Depth: Load

400 x 600mm: 5 KW

500 x 900mm: 8 KW

550 x 750mm: 8 KW

400 x 1100mm: 10 KW

500 x 1200mm: 10 KW

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003