

Product Name :
Horizontal Steam Autoclave

Product Code :
NLAB-LABORATORYEM47009



Description :

Horizontal Steam Autoclave

Technical Specification :

Horizontal Steam Autoclave -

Features:-

Pressure switch control.

Semi automatic operation or microprocessor based fully automatic operation.

High pressure high vacuum option.

Low water level.

Specifications:-

250 Ltrs

600x500x200mm
Dimension

92Kwh
Consumption


0.22 tkg/cm2
pressure

0.10 m3
Temperature

400 kg

800 kg
Weight

We are leading manufacturers, suppliers of Horizontal Steam Autoclave for Autoclave. Contact us to get high quality Horizontal Steam Autoclave for Autoclave for schools, colleges, universities, research labs, laboratories and various industries.



Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003