

Product Name :
Hospital Autopsy Anthropometric Rods

Product Code :
NLABMEDICAL-M2260001



Description :

Hospital Autopsy Anthropometric Rods

Technical Specification :

Anthropometric Rods are made of Brass.

These are used for vertical and horizontal measurements of long bones at the same time consisting of a set of four long and two small rods.

Four rods are attached at exclusive points to become a single rod to measure up to 200cm lengthwise.

Two small rods facilitate to clamp the bone and measure width-wise at the same time up to 25cm.

Chromium-plated brass rods have a smooth finish, movement, and clear markings.

We are Manufacturers, Exporters, and Suppliers of Anthropometric Rod at tender and export competitive prices from India.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com