

**Product Name :**  
Hospital Autopsy Anthropometric Rods

**Product Code :**  
NLABMEDICAL-M2260001



**Description :**

Hospital Autopsy Anthropometric Rods

**Technical Specification :**

Anthropometric Rods are made of Brass.

These are used for vertical and horizontal measurements of long bones at the same time consisting of a set of four long and two small rods.

Four rods are attached at exclusive points to become a single rod to measure up to 200cm lengthwise.

Two small rods facilitate to clamp the bone and measure width-wise at the same time up to 25cm.

Chromium-plated brass rods have a smooth finish, movement, and clear markings.

We are Manufacturers, Exporters, and Suppliers of Anthropometric Rod at tender and export competitive prices from India.

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)