

Product Name :
Hospital Electric Blood Donor Chair / Dialysis Chair

Product Code :
NLABMEDICAL-M2700014



Description :

Hospital Electric Blood Donor Chair / Dialysis Chair

Technical Specification :

Specification :-

Electric Blood chair for profiling of height, tilt, backrest and foot rest angles.

Three section table.

Adjustable armrests, footrest and backrests.

Height range of approx. 55cm-110cm.

Capacity of approx. 250 kg.

IV pole.

Dimensions approx : 210 L x 70 W x 55 H cms.

Optional Accessories at extra cost:-

Battery back-up.

Patient lap tray.

We are Manufacturer, Exporters and Suppliers of Electric Blood Donor Chair / Dialysis Chair at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003