

Product Name :
Hospital Full Fowler Bed ABS Top

Product Code :
NLABMEDICAL-M2660009



Description :

Hospital Full Fowler Bed ABS Top

Technical Specification :

Frame work of bed is made of Rectangular Tube.

4 Section top made of CRCA sheet uniformly perforated with latest parallel beam type surface.

Backrest and knee rest positions obtained by Smooth Crank System with foldable handles.

ABS Moulded Head and Foot Boards .

Bed is mounted on 100mm Dia Castor, 2 with brakes.

SS Telescopic I.V Rod with four location.

Holder for Urine Bag.

Product Size : 200 L x 90 W x 56 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Full Fowler Bed ABS Top at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003