

Product Name :
Hospital Full Fowler Bed SS Panel

Product Code :
NLABMEDICAL-M2660006



Description :

Hospital Full Fowler Bed SS Panel

Technical Specification :

Frame work of bed is made of Rectangular Tube.

4 Section top made of CRCA perforated Sheet.

Laminated SS Head and Foot panels.

Bed is mounted on 100mm Dia Castor, 2 with brakes.

Holder for Urine bag.

Backrest & Kneerest elevation positions obtained by Smooth Crank System with foldable handles.

4 Location for IV Rod.

Product Size: 200 L x 90 W x 56 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Full Fowler Bed SS Panel at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003