

**Product Name :**  
Hospital Full Fowler Bed Standard

**Product Code :**  
NLABMEDICAL-M2660005



**Description :**

Hospital Full Fowler Bed Standard

**Technical Specification :**

Frame work of bed is made of Rectangular Tube.

4 Section top made of CRCA perforated Sheet.

Tubular Head & Foot Bows are of unequal height with vertical tubular supports.

Bed is mounted on 100mm Dia Castor, 2 with brakes.

Holder for Urine bag.

Backrest & Kneerest elevation positions obtained by Smooth Crank System with foldable handles.

4 location for IV Rods.

Product Size: 200 L x 90 W x 56 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Full Fowler Bed Standard at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003