

<b>Product Name :</b> Hospital Instrument Cupboards	<b>Product Code :</b> NLABMEDICAL-M3270003
<b>Description :</b>	
Hospital Instrument Cupboards	
<b>Technical Specification :</b>	
Size: 960 x 400 x 1750 mm.	
Epoxy coated mild steel framework.	
Stainless steel bottom.	
Hospital Instrument Cabinet Hospital Instrument Cabinet are Size: 960 x 400 x 1750 mm	
Hospital Instrument Cabinet are Epoxy coated mild steel framework	
Hospital Instrument Cabinet are Stainless steel bottom	
Hospital Instrument Cabinet, Medical Cabinets Cupboards exporters and Medical Cabinets Cupboards manufacturers.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003