

Product Name :
Hospital Pediatric Bed Semi Fowler

Product Code :
NLABMEDICAL-M2660031



Description :

Hospital Pediatric Bed Semi Fowler

Technical Specification :

Frame is made of Rectangular pipe / Round Tube.

Top is made of perforated sheet.

It has 2 Section Top made of perforated Sheet.

Back Rest Section is manoeuvred by Screw Handle from Foot End.

Tubular Head & Foot Bows are of equal height with vertical tubular supports mounted on wheels.

The Bed is mounted on 10cm Dia Castors , 2 with brakes.

Safety Side railings can easily be drop down.

Product Size : 137 L x 76 W x 60 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Pediatric Bed Semi Fowler at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003