

**Product Name :**  
Hospital Semi Fowler Bed ABS Top

**Product Code :**  
NLABMEDICAL-M2660017



**Description :**

Hospital Semi Fowler Bed ABS Top

**Technical Specification :**

Frame work of bed is made of Rectangular Tube.

2 section top made of CRCA sheet uniformly perforated with latest parallel beam type surface.

Backrest positions obtained by Smooth Crank System with foldable handles.

ABS moulded Head and Foot panels .

Bed mounted on 100mm dia. swivel castors, 2 with brakes.

Product Size: 198 L x 90 W x 56 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Semi Fowler Bed ABS Top at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003