

Product Name :
Hospital Semi Fowler Bed Wire Mesh Top without Wheels

Product Code :
NLABMEDICAL-M2660014



Description :

Hospital Semi Fowler Bed Wire Mesh Top without Wheels

Technical Specification :

Frame work of bed is made of Rectangular Tube.

2 Section Top made of Welded Wire Mesh.

Backrest positions obtained by Smooth Crank System with handle.

Tubular Head & Foot Bows are of unequal height.

Legs mounted on hard rubber.

Product Size: 198 L x 90 W x 56 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Semi Fowler Bed Wire Mesh Top without Wheels at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003