

**Product Name :**  
Hospital Semi Fowler Bed Wire Mesh Top without Wheels

**Product Code :**  
NLABMEDICAL-M2660014



**Description :**

Hospital Semi Fowler Bed Wire Mesh Top without Wheels

**Technical Specification :**

Frame work of bed is made of Rectangular Tube.

2 Section Top made of Welded Wire Mesh.

Backrest positions obtained by Smooth Crank System with handle.

Tubular Head & Foot Bows are of unequal height.

Legs mounted on hard rubber.

Product Size: 198 L x 90 W x 56 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Semi Fowler Bed Wire Mesh Top without Wheels at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003