

**Product Name :**  
Hospital Waiting Chair Metal 2 Seater

**Product Code :**  
NLABMEDICAL-M2710005



**Description :**

Hospital Waiting Chair Metal 2 Seater

**Technical Specification :**

Frame is made of Round & Rectangular Tube.

The body section is mounted on 2 legs Round Tube Stable Base.

The waiting chair has Metal Perforated Seat with Back Rest Support.

It has Arm Support on both the corners.

We are Manufacturer, Exporters and Suppliers of Waiting Chair Metal 2 Seater at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003