

Product Name :
Hospital Waiting Chair Metal 3 Seater

Product Code :
NLABMEDICAL-M2710006



Description :

Hospital Waiting Chair Metal 3 Seater

Technical Specification :

Frame is made of Round & Rectangular Tube.

The body section is mounted on 2 legs Round Tube Stable Base.

The waiting chair has Metal Perforated Seat with Back Rest Support.

It has Arm Support on both the corners.

We are Manufacturer, Exporters and Suppliers of Waiting Chair Metal 3 Seater at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003