

Product Name :
Hospital Waste Bin Metal Perforated

Product Code :
NLABMEDICAL-M2720010



Description :

Hospital Waste Bin Metal Perforated

Technical Specification :

Waste bins comes without lid and are made of Stainless Steel.

These dustbins are in Perforated Sheet.

These are light weight, strong and durable.

These are hygienic and easy to clean.

These bins are rust free and maintainance free.

We are Manufacturer, Exporters and Suppliers of Waste Bin Metal Perforated at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

