

Product Name :
Hot Digital Water Bath 2 Liter

Product Code :
NLAB-EDUSCHOOLEQ45022



Description :

Hot Digital Water Bath 2 Liter

Technical Specification :

The digital water bath is a great economical solution for general purpose use at both high school and college level.

This is 2L with a temperature range of 5 to 99°C.

Capacity:

Covering

Voltage	
50-60Hz	
Consuming Power:	
400W	
Temperature Fluctuation:	
±0.1°C	
Temperature Range:	
5-99°C	
Chamber Dimensions:	
600x800x1000mm	
Exterior Dimensions:	
1515x915x1065mm	
Gross Weight:	
80kg	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003