

Product Name :
Hot Digital Water Bath 4 Liter

Product Code :
NLAB-EDUSCHOOLEQ45023



Description :

Hot Digital Water Bath 4 Liter

Technical Specification :

The digital water bath is a great economic solution for general purpose use at both high school and college level.
This is 4L with a temperature range of 5 to 99°C.

Capacity:

Coverings

Voltage	50-60Hz
Consuming Power:	600W
Temperature Fluctuation:	10-20°C
Temperature Range:	5-90°C
Internal Dimensions:	12x16x18"
External Dimensions:	20x20x10.5"
Gross Weight:	62kg±0.5kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003