

Product Name :
Hot Digital Water Bath 8 Liter

Product Code :
NLAB-EDUSCHOOLEQ45024



Description :

Hot Digital Water Bath 8 Liter

Technical Specification :

The digital water bath is a great economic solution for general purpose use at both high school and college level. This is 8L with a temperature range of 5 to 99°C.

Capacity:	
Color:	
Dimensions:	
Weight:	

Voltage	50-60Hz
Consuming Power:	
Temperature Fluctuation:	
Temperature Range:	
Chamber Dimensions:	
External Dimensions:	
Weight:	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003