

Product Name :
Humidifier

Product Code :
NLAB-SCIENTIFICEQ12014



Description :

Humidifier

Technical Specification :

Volume capacity: 40L

Display: High Contrast Display (HCD)

Tdew: 20...75°C

Pressure: 0...20 bar

Relative humidity: 10...100 %

Temperature: Optional 20...300 °C

Flow (l/min): 0...100 , 0...250, 0...500 ,0...1000

Visual alarm Indicator: Present

Thermal cycling system: Peltier-based, 96-well block

Sample Ramp Rate: +/-1.1°C/sec

Peak Block Ramp Rate: 1.6°C/sec

Temperature range: 4°C – 100°C

Temperature Accuracy: $\pm 0.25^\circ\text{C}$ of setpoint/display

temperature measured at 3 minutes after clock start

Temperature Uniformity: $\pm 0.50^\circ\text{C}$, 30 seconds after clock start

Power Connection: 240 V, 50/60 Hz;

Provide the warranty period of 24 months

Regulating humidity levels in the laboratory.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003