

Product Name :
Hydraulic Bench

Product Code :
NLAB-EDUCATIONTVET14001



Description :

Hydraulic Bench

Technical Specification :

Basic module for mounting fluid mechanics experiments (to be ordered separately)

Closed water circuit with submersible pump, storage tank and measuring tank.

The bench is constructed from lightweight reinforced plastic. The assembly is corrosion resistant and is mounted on casters.

An open channel is integrated into the bench table.

A set for volumetric measurements is included in the basic supply.

Technical characteristics:

Centrifugal pump :

Flow: max. 1.35 l / s

Pressure: max. 21 m of water

Power: 370 W

Storage tank capacity: 250 liters

High flow tank: 40 liters

Low flow tank: 6 liters

Working height: 1 meter

Dimensions (HxWxD): 1 x 1.13 x 0.73 m

Gross weight: 160 kg

Power supply: 220V - 50Hz

Recommended accessories:

Demo stopwatch, diameter 13 cm

Caliper in stainless steel 0 to 160 mm, precision 1/20 th

Wypall® roll in recyclable cellulose, 1 roll of 380 meters

Dispenser for paper roll up to 34 cm

Hose reel cart with 25m of garden hose and set of quick couplings

1/2 "tap with quick coupling

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003